

Medentika<sup>®</sup> Prosthetic Solutions

# NOVALOC<sup>®</sup> MEDENTILOC<sup>™</sup> DIGITAL WORKFLOW



# MEDENTIKA® NOVALOC® MATRIX SYSTEM

Novaloc® redefines denture retention with its precision-engineered attachment system. Designed with patient comfort in mind, the Novaloc® system facilitates convenient insertion and removal of prostheses, while allowing for optimal adjustment of denture retention according to different clinical needs.



## Divergence Compensation

In combination with the angled abutments, you can compensate for divergences of up to 70° between the implants.

## Gingival Heights

Both the straight and angled Novaloc® abutments are available in 5 different gingival heights.

## Angulation

15° angulation of the abutments.



## Matrix Housing

The matrix housing is available in titanium and also PEEK.



## Screw Head Opening

The small screw head opening of the straight Novaloc® abutment reduces food packing.



## ADLC Surface

The surface quality of the ADLC coating (amorphous diamond-like carbon) provides smooth surface and high hardness for wear resistance and low maintenance.



## Retention Insert

Made from PEEK high-performance plastic and manufactured with high precision, it can optimally absorb lateral pressure thanks to the patented design. Available in 6 different retention forces.

# MEDENTILOC™ ABUTMENTS



The MedentiLOC™ abutment offers high precision and is fully compatible with the Novaloc® matrix system. It provides a particularly cost-effective solution for overdenture retention. A notable advantage of the MedentiLOC™ abutment is its straightforward design, which allows for chairside treatment.




THE MEDENTILOC™ ABUTMENTS ARE AVAILABLE IN STRAIGHT AND ANGLED DESIGN AS WELL AS IN FIVE DIFFERENT GINGIVAL HEIGHTS.

## MULTI-UNIT SCREW PATRIX MEDENTILOC™ / NOVALOC®

The multi-unit screw patrices, when used with our multi-unit abutment, provide an additional option for restoration. They facilitate a seamless transition from a screw-retained restoration to a removable overdenture construction.

The multi-unit screw patrices, along with MedentiLOC™ and Novaloc® angled options, accommodate implant angulation of up to 30°.

These screw patrices are secured using the appropriate screwdriver for each implant system.



# NOVALOC®/MEDENTILOC™ DIGITAL

## DIGITAL PRECISION FOR REMOVABLE PROSTHETICS

Novaloc®/MedentiLOC™ Digital represents our innovation that advances removable prosthetics into the digital era of dentistry.

This fully digital workflow includes Scanbodies for the Novaloc® and MedentiLOC™ Abutments, sterilizable titanium Repositional Analog, and Retention Screws all designed for precision and efficiency.

With the library available for 3Shape and exocad, it simplifies data storage and sharing, eliminates conventional impression materials, and reduces patient visits. Novaloc®/MedentiLOC™ Digital is a versatile, indispensable tool for modern dental practices and is compatible with major Intra-oral Scanning systems and designed for ease of scanning.



### ABUTMENT LEVEL WORKFLOW




### ABUTMENT LEVEL WORKFLOW WITH PRINTED MODEL





# NOVALOC® / MEDENTILOC™ DIGITAL

## PRODUCT DESCRIPTION AND OVERVIEW



**Abutment Level Scanbody**

PEEK material for efficient scanning

**Repositional Analogs**

Represent the precise implant location

**Retention Screw for 3D Models**

Placement of the analog in the 3D-printed model is possible using click-fit" and screw retained method

# NOVALOC® / MEDENTILOC™ DIGITAL WORKFLOW

The virtually 100% digital restorative treatment workflow  
for implant-based removable prosthetic solution.



DATA ACQUISITION



TREATMENT PLANNING



IMPLANT SURGERY



ABUTMENT PLACEMENT

Novaloc® / MedentiLOC™



INTRAORAL SCAN

Abutment level scanbodies



PROSTHESIS DESIGN

Precise post-surgical design



WITHOUT MODEL


**Direct technique**

Intraoral abutment selection  
and workflow

WITH MODEL

**Direct or indirect technique**

Patient model creation with  
abutment analogs



LIBRARY AVAILABLE FOR  
3SHAPE AND EXOCAD

# BENEFITS OF DIGITAL WORKFLOW



## EFFICIENCY

- Precision
- Speed
- Chairside digital workflow
- Abutment library available for 3Shape and exocad



## FLEXIBILITY

- Works with major Intra-oral Scanning systems
- Easy to scan due to design and materials used



## PATIENT EXPERIENCE

- No conventional impression material is needed
- Less patient visits




## SCANBODY AND REPOSITIONAL ANALOG PRODUCT INFORMATION

ARTICLE NO'S	PRODUCT DESCRIPTION	MATERIAL
2010.726-NOV	Novaloc® Scanbody 4 pack	PEEK
2010.727-NOV	Novaloc® Scanbody 10 pack	PEEK
2010.728-NOV	Novaloc® Model Analog CAD/CAM 4 pack	Titanium
2010.729-NOV	Novaloc® Model Analog CAD/CAM 10 pack	Titanium

SCAN OR CLICK FOR OUR ONLINE SHOP:




SCAN OR CLICK TO DOWNLOAD OUR LIBRARIES:





# STEP-BY-STEP APPLICATION



1




- Ensure that the seating surface of the abutment is clean, free of residue, and dry.
- Place the Scanbody on each Novaloc® / MedentiLOC™ abutment and press down firmly. (If there is a gap between the scanbody and the gingiva, Novaloc® / MedentiLOC™ mounting collar can be placed on the abutments to close the gap for better scanning performance)

2



- Continue scanning by following the instructions of the software and hardware IFU.
- Select “just scan” or denture workflow in the scanning software.

3




- Process the scan and send it to the dental laboratory for the design and manufacturing of the final prosthesis.
- The Scanbody creates a recess in the prosthesis base to allow the housing to be seated without interference from the surrounding prosthesis acrylic.

SCAN OR CLICK TO SEE  
THE STEP BY STEP VIDEO  
IN **3SHAPE**



SCAN OR CLICK TO SEE  
THE STEP BY STEP VIDEO  
IN **EXOCAD**






# BEYOND THE ORIGINAL



CLICK OR SCAN THE CODE FOR  
MORE DETAILS ABOUT OUR  
PRODUCTS IN OUR ESHOP





# RE-ENGINEERED COMPATIBLE PROSTHETICS.



A Straumann Group Brand

© MEDENTiKA® GmbH, 2024. All rights reserved.  
MEDENTiKA® and/or other trademarks and logos from MEDENTiKA® mentioned herein are the trademarks or registered trademarks of MEDENTiKA® GmbH and/or its affiliates.

Please note that this brochure contains images created by AI.

PM02\_01\_0064\_EN\_09/24

**MEDENTiKA® GmbH**  
Hammweg 8–10  
76549 Hügelsheim  
Germany  
info@medentika.de  
www.medentika.com

